

Exclaimer Mail Archiver 3.01

Release Notes

Release Number

The current Release Number for this product is: 3.01.60130.0

System Requirements

Hardware	Minimum Requirements	Recommended Requirements
CPU	x64 architecture-based computer. Intel Itanium family IA64 processor are not supported.	x64 architecture-based computer. Intel Itanium family IA64 processor are not supported.
Memory	4GB	8GB
Disk Space for program files	350 MB	350MB
Screen Resolution	1024 x 768 pixels	1280 x 1024 or higher
Software	Supported Versions	
Operating Systems	Windows Server 2008 x64 Windows Server 2008 SP1 x64 Windows Server 2008 SP2 x64 Windows Server 2008 R2 x64 Windows Small Business Server 2008 Windows Small Business Server 2011 Windows Vista x64 Windows Vista SP1 x64 Windows Vista SP2 x64 Windows 7 x64 Windows 7 SP1 x64 Windows 8 x64 Windows Server 2012 (Recommended)	
Exchange Server	Microsoft Exchange Server 2007 SP1 Rollup 5 or higher Microsoft Exchange Server 2007 SP2 Microsoft Exchange Server 2007 SP3 Microsoft Exchange Server 2010 Microsoft Exchange Server 2010 SP1 Microsoft Exchange Server 2010 SP2 Microsoft Exchange Server 2013*	
Microsoft .Net Framework	Microsoft .Net Framework 3.5 SP1 and above.	

*Exclaimer Mail Archiver should not be installed directly on an Exchange 2013 server.

The following software is required for Exclaimer Mail Archiver. These components will be installed via the Setup program if not present on the system at the time of installation.

Software (Prerequisites)	Minimum Requirements
Microsoft IIS	Internet Information Services 7 and above
Microsoft ASP.NET MVC	Microsoft ASP.NET MVC 2

Note: Desktop operating systems (Windows Vista, Windows 7 and Windows 8) have a maximum of 10 concurrent connections for IIS. Exclaimer recommends using these operating systems for evaluation purposes only.

Downloading

Please refer to the following link on the Exclaimer website:

<http://www.exclaimer.com/mail-archiver/download>

Installing

Download the 'Mail Archiver 3.0' installation file 'setup.exe' and run on one of the supported operating systems listed under the Software heading in the System Requirements section. An MSI is available from support by contacting support@exclaimer.com.

Uninstalling

Uninstall can be achieved via the 'setup.exe' or MSI for the currently installed version of 'Mail Archiver 3.0' or via 'Programs and Features' in your operating system.

The configuration files are located in "\\ProgramData\\Exclaimer Ltd\\Mail Archiver" and will not be removed as part of the uninstall process.

The uninstall process will not lose any of your configuration or settings. However as a backup precaution you may wish to export your configuration.

Export Steps

1. Open the Exclaimer Console.
2. Select the 'Exclaimer' node.
3. Open the 'Action' menu.
4. Select 'Export Configuration...'
5. Select a location and filename for the export.
6. Press Save.

Note: This will export all your settings and licensing data only. Your archived data will not be exported.

Upgrading

Upgrading can be achieved by following the steps above to uninstall the current version and then install the new version. Journal rules will not be disabled, messages will still be delivered to the Journal Mailbox. Journal Mailbox Importers will continue to archive messages once the installation of the new version has been completed.

Fixed Issues

- "Timezone ID was not found" error when connecting to the web search sign-in page when Mail Archiver is installed on certain non-English versions of Windows.
- Several display issues with column sorting when the search query cache has expired.
- Journal Mailbox Importer cannot connect when Journal Inbox name is not "Inbox".
- Unable to import configuration settings when auditing feature is enabled.